

Risk Management 2.0



2017 annual INCOSE international workshop

Los Angeles, CA, USA January 28 - 31, 2017

MBSE - Transforming Medical Device Design

Casey Medina & Kristina Fuerst



5 - 1-1FH

Casey Medina, CSEP

Casey currently practices systems engineering in the medical device industry. He has developed products in the areas of automated stem cell growth, patient monitoring, pulse-oximetry, blood collections and therapeutics. Casey holds multiple patents in the medical device sector. Professionally, he is focused on developing the art of systems engineering in a manner that fosters adoption and acceptance by organizations resistant to change. He currently applies systems engineering practices and principles to enable rapid medical device development and process design and is working to enhance the use of MBSE as an enabler for usability and human factors analyses.







Kristina Fuerst, ASEP

Kristina has been working in the medical device industry as a systems engineer since 2012 and is responsible for the requirements engineering and system architecture development for her team. Kristina developed and implemented an intensive, two-day requirements training course with Casey. Kristina graduated from the University of Colorado at Boulder with a B.S. in Chemical and Biological Engineering and is pursuing her Master's degree in Biomedical Engineering with Rutgers University. She received her ASEP in 2013.


















Terumo BCT

Established:	1964 as COBE Laboratories	
President and CEO:	David Perez	
Number of Associates:	5,300+	
Global Headquarters:	Lakewood, Colorado	
Regional Headquarters:	Brussels, Buenos Aires, Hong Kong, Tokyo	
Manufacturing Facilities:	United States, Belgium, Japan, Northern Ireland, India, Vietnam	
Main Business:	A global leader in blood component, therapeutic apheresis and cellular technologies, and the only company with the unique combination of apheresis collections, manual and automated whole blood processing, and pathogen reduction. We believe in the potential of blood to do even more for patients than it does today. This belief inspires our innovation and strengthens our collaboration with customers.	
Parent Company:	Terumo Corporation; Tokyo Stock Exchange (4543)	

Terumo BCT





SysML Elements and Diagrams Created with MagicDraw by No Magic

Proceedings of the 25th Annual International Symposium


"From Asking Forgiveness to Saying You're Welcome, Introducing Requirements Engineering to Medical Device Development" Casey Medina, Kristina Fuerst

"The "V" Model Reloaded" Casey Medina 33%

Where We Started

- No understanding of system architecture outside of limited, localized, software practices
- No process for developing system architecture
- Design information largely captured as "tribal knowledge"

MBSE Coordinates Development Processes



Our Challenge

- "This is telling me how to do my job!"
- "We know this already."
- "This isn't possible to know."
- "This is too much information."
- Perception that schedule pressure won't allow this work to happen
- Definition of Intended use (mission statement)
- Process for determining Objectives and Goals
- Process for creating User/Env analysis
 - interactions
 - characteristics
- Descriptions of Situational Behaviors
- Definitions of System Behavior
- -TBCT Associate

- Teach system modeling to encourage its integration in development processes
 - Examples of various system operational scenarios
 - Actual product example

Fnablers

- SAD template
- Example description of system behavior
- Architecture guidance, publications, best practices, etc.

- 1 Demonstrate ability to generate at least 2 levels of hierarchy
- 2 Understand what is meant by system behavior and how it is different from system performance
- 3 Understand hierarchy and flow of system behaviors
- 4 Differentiate black and white box behaviors
- 5 Demonstrate ability to create FFBD and context diagrams for system behavior from intended use
- 6 Understand how to use the syntax to document the models
- 7 Understand these models are the backbone of the project. Dictates Req's, V&V risk, all parts of design

Models Answer Questions

Answers are expensive. Good answers cost the same. - Rick Griffith Teams are more comfortable solving problems than defining them.

Asking the questions right is as important as asking the right questions.

- Ask questions in a personal, engaging manner
- Listen intently be present
- Create a safe environment
 Help your team capitalize
- Help your team capitalize on opportunities to increase their knowledge of their system

How reall down this result in purior

How does the system respond if that happens?

Do I undustand all the interfaces?

What do the sales think of this design.

Did I capture the requirements accountally?

Here much coverage do I have on my VAV unitvities?

What do I roud hope to loan

What tighe of shaded will be 12 help the lease?

How many levels of Hancodry and necessary

How will I varify and solidate my model?

How self I share the results of my trodeling activities

Answers are expensive. Good answers cost the same.

- Rick Griffith

Teams are more comfortable solving problems than defining them.

Asking the questions right is as important as asking the right questions.

- Ask questions in a personal, engaging manner
- Listen intently be present
- Create a safe environment
- Help your team capitalize on opportunities to increase their knowledge of their system

How well does this need to perform?

How much should it weigh?

How does the system respond if that happens?

Do I understand all the interfaces?

What do the users think of this design?

Did I capture the requirements accurately?

How much coverage do I have on my V&V activities?

What do I need/hope to learn?

What type of model will best help me learn?

How much information do I need to represent?

How many levels of hierarchy are necessary?

How will I verify and validate my model?


How will I share the results of my modeling activities?

Structure Implements Function


































Segrit file is being debetom ung genetien betreffen in der de			
Securities to the continue of	rakondor	Agentowiti	ditain jugges meets never per neutram (p.1)
The control of the co	trafficers or terior	Noon Cast to be refinite one at rable. Noon Cast to common emission abstrage space. Noon Cast to complex as independent beyond.	
Amen fall also also also also also also also a	Medical Protocoled		- News phone waters to be as a side
Amon Africa de la como antica de mano apparen la como apparen	Minde		
Amount of the fail for the management of the property of the p	Name Translationer		
The second secon	Special Contract Con-		
And the second s	tru.		
-Moots monitorise be visually appoints and instillMoots to selable to use the IQM. The file own se		macrostic . Associated and the state of the	annother memoryoparatic regions for given and memory and properties of the second of t
		Roots monitor selbs vassells apposing and instill	- floors to be able to use the Cife. Theolike own of


- Duty cycle
- "rise and fall" times dictated by centrifuge speed
- Necessary life of LED
- Camera sensitivity
- Acceptable light levels for camera




- Centrifuge position
- LED expected life











Casey's Oversimplified-but-Really-Useful Interface Model





Before MBSE approach, one product family contained ~30 "requirements" scattered about its documentation

After MBSE approach, that same product family contained ~ 350 requirements, ordered hierarchically.

The team identified V&V gaps and worked hard to fill them.

The team has a great baseline from which to make future enhancements or to deal with obselescence.

Making a "Quantum" Leap





















Controls

- Regulatory
 - -GMP requirements (closed system, traceability, etc.)
 - -Medical device standards and guidelines
- Hollow fiber technology constraint adherence of cells
- · Society's familiarity with flask culture
- · Physiological needs to cells
 - -pH
 - -temperature
- Contamination risks
 - -nonmutagenic (leachables, extractables)
 - -efficacy of cells
- Ability to isolate desired cell types from mixed populations
- Targeted cell types
- Ability of operator to create customized protocols/parameters
- · Automation constraints
- Customer locations (APAC, Europe, US, etc.)
 - -Differences in interface types (gas, power, etc.)
 - Acquisition needs (CE marking, etc.)
- Ability of system to detect anomalous conditions

Activity

Grow adherent and suspension cells in laboratory, clinical, and manufacturing environments

Inputs

- Initial population of cells in solution (isolated or mixed population)
- Cell growth nutrients/fluids
- · Cell growth system
- · Unique identification of expansion cycle
- Oxygen supply (gas supply)
- Unique identification for cell expansion
- User interaction data

Enablers

- Skilled personnel (operators, maintainers, trainers)
- System maintenance
- Sterile culture environment (single use set)
- · Closed culture environment
- Trained operator
- Cell metabolic rate monitoring method
- Automated protocols
- Data information storage system
- User instructions (service, installation, operation)
- Power supply
- Gas control mechanism
- Cellular respiration byproduct waste management method

Outputs

· Final population of cells


Cell number


Viability


Potency


Cellular products (proteins, vectors, etc.)


- Waste materials
 - -Fluid waste
 - -Cell respiratory waste (CO2, lactate, etc.)
 - -Spent single-use set
 - -Dead/unviable cells
- · Exhausted cell growth system
- Reports
 - -User Interaction and system state linked to unique procedure identification
 - -Troubleshooting information (servicing purposes)
- · Alerts of anomalous conditions





























Risk Management 2.0


























Guiding Thoughts

- Keep it simple to give teams confidence and practice in applying fundamentals. Don't overwhelm them!
- Modeling techniques and methodologies may not always be applied "by the book."
- Help teams understand why they're building the model to help them focus and drive for results.
- Teams need to feel that the model helps them get to the finish line more efficiently.
- Keep the end in mind.

Please Remember:

- What is obvious to an SE is not necessarily obvious to others
- Empower your teams to embrace the process
- Allow your teams time to digest and understand the material

We're only successful if we enable our teams to be successful!